

Antibodies Product Datasheet

Product Name:	EPHA5 (Phospho-Tyr650) Antibody
Catalog Number:	A9829
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHA5 (Phospho-Tyr650)
Concentration:	1 mg/ml
Gene Name:	EPHA5 BSK EHK1 HEK7 TYRO4
Protein Name:	EPHA5 (Phospho-Tyr650)
Other Name:	Ephrin type-A receptor 5 (EC 2.7.10.1) (Brain-specific kinase) (EPH homology kinase 1) (EHK-1) (EPH-like kinase 7) (EK7) (hEK7)
Human Swiss Prot No.:	P54756
Formulation:	Liquid in PBS containing 50% glycerol, and 0.174% sodium azide.
Specificity:	This antibody detects endogenous levels of EPHA5 (Phospho-Tyr650) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year