

Antibodies Product Datasheet

Product Name:	EPHB1/2/3/4 (Phospho-Tyr600/602/614/596) Antibody
Catalog Number:	A9833
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHB1/2/3/4 (Phospho-Tyr600/602/614/596)
Concentration:	1 mg/ml
Gene Name:	EPHB1 ELK EPHT2 HEK6 NET
Protein Name:	EPHB1/2/3/4 (Phospho-Tyr600/602/614/596)
Other Name:	Ephrin type-B receptor 1 (EC 2.7.10.1) (ELK) (EPH tyrosine kinase 2) (EPH-like kinase 6) (EK6) (hEK6) (Neuronally-expressed EPH-related tyrosine kinase) (NET) (Tyrosine-protein kinase receptor EPH-2)
Human Swiss Prot No.:	P54762
Formulation:	Liquid in PBS containing 50% glycerol, and 0.187% sodium azide.
Specificity:	This antibody detects endogenous levels of EPHB1/2/3/4 (Phospho-Tyr600/602/614/596) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year