

Antibodies Product Datasheet

Product Name:	Gz- α (Phospho-Ser16) Antibody
Catalog Number:	A9835
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Gz- α (Phospho-Ser16)
Concentration:	1 mg/ml
Gene Name:	GNAZ
Protein Name:	Gz- α (Phospho-Ser16)
Other Name:	Guanine nucleotide-binding protein G(z) subunit alpha (G(x) alpha chain) (Gz-alpha)
Human Swiss Prot No.:	P19086
Formulation:	Liquid in PBS containing 50% glycerol, and 0.192% sodium azide.
Specificity:	This antibody detects endogenous levels of Gz- α (Phospho-Ser16) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year