

Antibodies Product Datasheet

Product Name:	LIMK1 (Phospho-Ser310) Antibody
Catalog Number:	A9837
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LIMK1 (Phospho-Ser310)
Concentration:	1 mg/ml
Gene Name:	LIMK1 LIMK
Protein Name:	LIMK1 (Phospho-Ser310)
Other Name:	LIM domain kinase 1 (LIMK-1) (EC 2.7.11.1)
Human Swiss Prot No.:	P53667
Formulation:	Liquid in PBS containing 50% glycerol, and 0.204% sodium azide.
Specificity:	This antibody detects endogenous levels of LIMK1 (Phospho-Ser310) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year