

Antibodies Product Datasheet

Product Name:	lipin1 Phospho-Ser889 Antibody
Catalog Number:	A9841
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human lipin1 ser889
Concentration:	1 mg/ml
Gene Name:	LPIN1 KIAA0188
Protein Name:	lipin1 ser889
Other Name:	Phosphatidate phosphatase LPIN1 (EC 3.1.3.4) (Lipin-1)
Human Swiss Prot No.:	Q14693
Formulation:	Liquid in PBS containing 50% glycerol, and 0.223% sodium azide.
Specificity:	This antibody detects endogenous levels of lipin1 ser889 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year