

Antibodies Product Datasheet

Product Name:	MEF-2 Antibody
Catalog Number:	B0019
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEF2A. AA range:374-423
Concentration:	1 mg/ml
Gene Name:	MEF2A
Protein Name:	Myocyte-specific enhancer factor 2A
Other Name:	MEF2A; MEF2; Myocyte-specific enhancer factor 2A; Serum response factor-like protein 1
Human Swiss Prot No.:	Q02078
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEF-2 Polyclonal Antibody detects endogenous levels of MEF-2 protein.
Storage and Stability:	-20°C/1 year