

Antibodies Product Datasheet

Product Name:	c-Src Antibody
Catalog Number:	B0045
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from c-Src . at AA range: 360-440
Concentration:	1 mg/ml
Gene Name:	SRC
Protein Name:	Proto-oncogene tyrosine-protein kinase Src
Other Name:	SRC; SRC1; Proto-oncogene tyrosine-protein kinase Src; Proto-oncogene c-Src; pp60c-src; p60-Src
Human Swiss Prot No.:	P12931
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	c-Src Polyclonal Antibody detects endogenous levels of c-Src protein.
Storage and Stability:	-20°C/1 year