

Antibodies Product Datasheet

Product Name:	LIMK-1/2 Antibody
Catalog Number:	B0047
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LIMK1/2. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	LIMK1/LIMK2
Protein Name:	LIM domain kinase 1/LIM domain kinase 2
Other Name:	LIMK1; LIMK; LIM domain kinase 1; LIMK-1; LIMK2; LIM domain kinase 2; LIMK-2
Human Swiss Prot No.:	P53667/P53671
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LIMK-1/2 Polyclonal Antibody detects endogenous levels of LIMK-1/2 protein.
Storage and Stability:	-20°C/1 year