

Antibodies Product Datasheet

Product Name:	GRF-1 Antibody
Catalog Number:	B0433
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRF-1. AA range:1071-1120
Concentration:	1 mg/ml
Gene Name:	ARHGAP35
Protein Name:	Rho GTPase-activating protein 35
Other Name:	ARHGAP35; GRF1; GRLF1; KIAA1722; Rho GTPase-activating protein 35; Glucocorticoid receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A
Human Swiss Prot No.:	Q9NRY4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRF-1 Polyclonal Antibody detects endogenous levels of GRF-1 protein.
Storage and Stability:	-20°C/1 year