

Antibodies Product Datasheet

Product Name:	KOR-1 Antibody
Catalog Number:	B0448
Reactivity:	Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from mouse KOR-1. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	OPRK1
Protein Name:	Kappa-type opioid receptor
Other Name:	OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KOR-1 Polyclonal Antibody detects endogenous levels of KOR-1 protein.
Storage and Stability:	-20°C/1 year