

## Antibodies Product Datasheet

|                               |                                                                                                                                                                           |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | NMDAε2 Antibody                                                                                                                                                           |
| <b>Catalog Number:</b>        | B0526                                                                                                                                                                     |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                         |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                                                    |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/5000. Not yet tested in other applications.                                                |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                          |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                         |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                     |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human NMDAR2B. AA range:1435-1484                                                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                   |
| <b>Gene Name:</b>             | GRIN2B                                                                                                                                                                    |
| <b>Protein Name:</b>          | Glutamate [NMDA] receptor subunit epsilon-2                                                                                                                               |
| <b>Other Name:</b>            | GRIN2B; NMDAR2B; Glutamate [NMDA] receptor subunit epsilon-2; N-methyl D-aspartate receptor subtype 2B; NMDAR2B; NR2B; N-methyl-D-aspartate receptor subunit 3; NR3; hNR3 |
| <b>Human Swiss Prot No.:</b>  | Q13224                                                                                                                                                                    |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                   |
| <b>Specificity:</b>           | NMDAε2 Polyclonal Antibody detects endogenous levels of NMDAε2 protein.                                                                                                   |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                              |