

Antibodies Product Datasheet

Product Name:	ACK Antibody
Catalog Number:	B0604
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACK1. AA range:250-299
Concentration:	1 mg/ml
Gene Name:	TNK2
Protein Name:	Activated CDC42 kinase 1
Other Name:	TNK2; ACK1; Activated CDC42 kinase 1; ACK-1; Tyrosine kinase non-receptor protein 2
Human Swiss Prot No.:	Q07912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACK Polyclonal Antibody detects endogenous levels of ACK protein.
Storage and Stability:	-20°C/1 year