

Antibodies Product Datasheet

Product Name:	eIF4G Antibody
Catalog Number:	B0642
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human eIF4G. AA range:1074-1123
Concentration:	1 mg/ml
Gene Name:	EIF4G1
Protein Name:	Eukaryotic translation initiation factor 4 gamma 1
Other Name:	EIF4G1; EIF4F; EIF4G; EIF4GI; Eukaryotic translation initiation factor 4 gamma 1; eIF-4-gamma 1; eIF-4G 1; eIF-4G1; p220
Human Swiss Prot No.:	Q04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF4G Polyclonal Antibody detects endogenous levels of eIF4G protein.
Storage and Stability:	-20°C/1 year