

Antibodies Product Datasheet

Product Name:	Ataxin-1 Antibody
Catalog Number:	B0771
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Ataxin 1. AA range:742-791
Concentration:	1 mg/ml
Gene Name:	ATXN1
Protein Name:	Ataxin-1
Other Name:	ATXN1; ATX1; SCA1; Ataxin-1; Spinocerebellar ataxia type 1 protein
Human Swiss Prot No.:	P54253
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ataxin-1 Polyclonal Antibody detects endogenous levels of Ataxin-1 protein.
Storage and Stability:	-20°C/1 year