

Antibodies Product Datasheet

Product Name:	CENP-A Antibody
Catalog Number:	B0858
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Centromeric Protein A. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CENPA
Protein Name:	Histone H3-like centromeric protein A
Other Name:	CENPA; Histone H3-like centromeric protein A; Centromere autoantigen A; Centromere protein A; CENP-A
Human Swiss Prot No.:	P49450
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CENP-A Polyclonal Antibody detects endogenous levels of CENP-A protein.
Storage and Stability:	-20°C/1 year