

Antibodies Product Datasheet

Product Name:	DAPK2 Antibody
Catalog Number:	B0899
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DAPK2. AA range:284-333
Concentration:	1 mg/ml
Gene Name:	DAPK2
Protein Name:	Death-associated protein kinase 2
Other Name:	DAPK2; Death-associated protein kinase 2; DAP kinase 2; DAP-kinase-related protein 1; DRP-1
Human Swiss Prot No.:	Q9UIK4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DAPK2 Polyclonal Antibody detects endogenous levels of DAPK2 protein.
Storage and Stability:	-20°C/1 year