

Antibodies Product Datasheet

Product Name:	BAM32 Antibody
Catalog Number:	B0901
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DAPP1. AA range:105-154
Concentration:	1 mg/ml
Gene Name:	DAPP1
Protein Name:	Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide
Other Name:	DAPP1; BAM32; HSPC066; Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide; hDAPP1; B lymphocyte adapter protein Bam32; B-cell adapter molecule of 32 kDa
Human Swiss Prot No.:	Q9UN19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAM32 Polyclonal Antibody detects endogenous levels of BAM32 protein.
Storage and Stability:	-20°C/1 year