

Antibodies Product Datasheet

Product Name:	IFN- α / β R α Antibody
Catalog Number:	B0952
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Interferon-alpha/beta Receptor alpha chain. AA range:436-485
Concentration:	1 mg/ml
Gene Name:	IFNAR1
Protein Name:	Interferon alpha/beta receptor 1
Other Name:	IFNAR1; IFNAR; Interferon alpha/beta receptor 1; IFN-R-1; IFN-alpha/beta receptor 1; Cytokine receptor class-II member 1; Cytokine receptor family 2 member 1; CRF2-1; Type I interferon receptor 1
Human Swiss Prot No.:	P17181
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFN- α / β R α Polyclonal Antibody detects endogenous levels of IFN- α / β R α protein.
Storage and Stability:	-20°C/1 year