

Antibodies Product Datasheet

Product Name:	Hic-5 Antibody
Catalog Number:	B1036
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Hic-5. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	TGFB1I1
Protein Name:	Transforming growth factor beta-1-induced transcript 1 protein
Other Name:	TGFB1I1; ARA55; Transforming growth factor beta-1-induced transcript 1 protein; Androgen receptor coactivator 55 kDa protein; Androgen receptor-associated protein of 55 kDa; Hydrogen peroxide-inducible clone 5 protein; Hic-5
Human Swiss Prot No.:	O43294
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hic-5 Polyclonal Antibody detects endogenous levels of Hic-5 protein.
Storage and Stability:	-20°C/1 year