

Antibodies Product Datasheet

Product Name:	HMG-14 Antibody
Catalog Number:	B1040
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HMG14. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	HMGN1
Protein Name:	Non-histone chromosomal protein HMG-14
Other Name:	HMGN1; HMG14; Non-histone chromosomal protein HMG-14; High mobility group nucleosome-binding domain-containing protein 1
Human Swiss Prot No.:	P05114
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HMG-14 Polyclonal Antibody detects endogenous levels of HMG-14 protein.
Storage and Stability:	-20°C/1 year