

Antibodies Product Datasheet

Product Name:	IL-7R Antibody
Catalog Number:	B1065
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-7R/CD127. AA range:410-459
Concentration:	1 mg/ml
Gene Name:	IL7R
Protein Name:	Interleukin-7 receptor subunit alpha
Other Name:	IL7R; Interleukin-7 receptor subunit alpha; IL-7 receptor subunit alpha; IL-7R subunit alpha; IL-7R-alpha; IL-7RA; CDw127; CD antigen CD127
Human Swiss Prot No.:	P16871
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-7R Polyclonal Antibody detects endogenous levels of IL-7R protein.
Storage and Stability:	-20°C/1 year