

Antibodies Product Datasheet

Product Name:	KCNQ2/3/4/5 Antibody
Catalog Number:	B1089
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Kv7.3/KCNQ3. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	KCNQ2
Protein Name:	Potassium voltage-gated channel subfamily KQT member 2
Other Name:	KCNQ2; Potassium voltage-gated channel subfamily KQT member 2; KQT-like 2; Neuroblastoma-specific potassium channel subunit alpha KvLQT2; Voltage-gated potassium channel subunit Kv7.2; KCNQ3; Potassium voltage-gated channel subfamily KQT me
Human Swiss Prot No.:	O43526/O43525/P56696/Q9NR82
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNQ2/3/4/5 Polyclonal Antibody detects endogenous levels of KCNQ2/3/4/5 protein.
Storage and Stability:	-20°C/1 year