

Antibodies Product Datasheet

Product Name:	MAP4 Antibody
Catalog Number:	B1092sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MAP4 Polyclonal
Concentration:	1 mg/ml
Gene Name:	MAP4
Protein Name:	Microtubule-associated protein 4 (MAP-4)
Other Name:	Microtubule-associated protein 4 (MAP-4)
Human Swiss Prot No.:	P27816
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAP4.
Storage and Stability:	-20°C/1 year