

Antibodies Product Datasheet

Product Name:	NIFK Antibody
Catalog Number:	B1106
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NIFK. AA range:200-249
Concentration:	1 mg/ml
Gene Name:	MKI67IP
Protein Name:	MKI67 FHA domain-interacting nucleolar phosphoprotein
Other Name:	MKI67IP; NIFK; NOPP34; MKI67 FHA domain-interacting nucleolar phosphoprotein; Nucleolar phosphoprotein Nopp34; Nucleolar protein interacting with the FHA domain of pKl-67; hNIFK
Human Swiss Prot No.:	Q9BYG3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NIFK Polyclonal Antibody detects endogenous levels of NIFK protein.
Storage and Stability:	-20°C/1 year