

Antibodies Product Datasheet

Product Name:	C3G Antibody
Catalog Number:	B1120
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAPGEF1. AA range:470-519
Concentration:	1 mg/ml
Gene Name:	RAPGEF1
Protein Name:	Rap guanine nucleotide exchange factor 1
Other Name:	RAPGEF1; GRF2; Rap guanine nucleotide exchange factor 1; CRK SH3-binding GNRP; Guanine nucleotide-releasing factor 2; Protein C3G
Human Swiss Prot No.:	Q13905
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C3G Polyclonal Antibody detects endogenous levels of C3G protein.
Storage and Stability:	-20°C/1 year