

Antibodies Product Datasheet

Product Name:	ROM-K Antibody
Catalog Number:	B1121
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ROMK/Kir1.1. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	KCNJ1
Protein Name:	ATP-sensitive inward rectifier potassium channel 1
Other Name:	KCNJ1; ROMK1; ATP-sensitive inward rectifier potassium channel 1; ATP-regulated potassium channel ROM-K; Inward rectifier K(+) channel Kir1.1; Potassium channel; inwardly rectifying subfamily J member 1
Human Swiss Prot No.:	P48048
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ROM-K Polyclonal Antibody detects endogenous levels of ROM-K protein.
Storage and Stability:	-20°C/1 year