

Antibodies Product Datasheet

Product Name:	FANCA Antibody
Catalog Number:	B1148
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FANCA. AA range:1121-1170
Concentration:	1 mg/ml
Gene Name:	FANCA
Protein Name:	Fanconi anemia group A protein
Other Name:	FANCA; FAA; FACA; FANCH; Fanconi anemia group A protein; Protein FACA
Human Swiss Prot No.:	O15360
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FANCA Polyclonal Antibody detects endogenous levels of FANCA protein.
Storage and Stability:	-20°C/1 year