

Antibodies Product Datasheet

Product Name:	IGF-IIR Antibody
Catalog Number:	B1211
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IGF2R. AA range:2381-2430
Concentration:	1 mg/ml
Gene Name:	IGF2R
Protein Name:	Cation-independent mannose-6-phosphate receptor
Other Name:	IGF2R; MPRI; Cation-independent mannose-6-phosphate receptor; CI Man-6-P receptor; CI-MPR; M6PR; 300 kDa mannose 6-phosphate receptor; MPR 300; Insulin-like growth factor 2 receptor; Insulin-like growth factor II receptor; IGF-II receptor;
Human Swiss Prot No.:	P11717
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGF-IIR Polyclonal Antibody detects endogenous levels of IGF-IIR protein.
Storage and Stability:	-20°C/1 year