

Antibodies Product Datasheet

Product Name:	Atm Antibody
Catalog Number:	B7018
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATM. AA range:1950-1999
Concentration:	1 mg/ml
Gene Name:	ATM
Protein Name:	Serine-protein kinase ATM
Other Name:	ATM; Serine-protein kinase ATM; Ataxia telangiectasia mutated; A-T mutated
Human Swiss Prot No.:	Q13315
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Atm Polyclonal Antibody detects endogenous levels of Atm protein.
Storage and Stability:	-20°C/1 year