

Antibodies Product Datasheet

Product Name:	Cdc2 Antibody
Catalog Number:	B7036
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC2. AA range:5-54
Concentration:	1 mg/ml
Gene Name:	CDK1
Protein Name:	Cyclin-dependent kinase 1
Other Name:	CDK1; CDC2; CDC28A; CDKN1; P34CDC2; Cyclin-dependent kinase 1; CDK1; Cell division control protein 2 homolog; Cell division protein kinase 1; p34 protein kinase
Human Swiss Prot No.:	P06493
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc2 Polyclonal Antibody detects endogenous levels of Cdc2 protein.
Storage and Stability:	-20°C/1 year