

Antibodies Product Datasheet

Product Name:	G3BP1 Antibody
Catalog Number:	B7089
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human G3BP-1. AA range:199-248
Concentration:	1 mg/ml
Gene Name:	G3BP1
Protein Name:	Ras GTPase-activating protein-binding protein 1
Other Name:	G3BP1; G3BP; Ras GTPase-activating protein-binding protein 1; G3BP-1; ATP-dependent DNA helicase VIII; hDH VIII; GAP SH3 domain-binding protein 1
Human Swiss Prot No.:	Q13283
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	G3BP1 Polyclonal Antibody detects endogenous levels of G3BP1 protein.
Storage and Stability:	-20°C/1 year