

Antibodies Product Datasheet

Product Name:	IGF-IR Antibody
Catalog Number:	B7115
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IGF1R. AA range:1131-1180
Concentration:	1 mg/ml
Gene Name:	IGF1R
Protein Name:	Insulin-like growth factor 1 receptor
Other Name:	IGF1R; Insulin-like growth factor 1 receptor; Insulin-like growth factor I receptor; IGF-I receptor; CD antigen CD221
Human Swiss Prot No.:	P08069
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGF-IR Polyclonal Antibody detects endogenous levels of IGF-IR protein.
Storage and Stability:	-20°C/1 year