

Antibodies Product Datasheet

Product Name:	LIMK-2 Antibody
Catalog Number:	B8045
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LIMK2. AA range:249-298
Concentration:	1 mg/ml
Gene Name:	LIMK2
Protein Name:	LIM domain kinase 2
Other Name:	LIMK2; LIM domain kinase 2; LIMK-2
Human Swiss Prot No.:	P53671
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LIMK-2 Polyclonal Antibody detects endogenous levels of LIMK-2 protein.
Storage and Stability:	-20°C/1 year