

Antibodies Product Datasheet

Product Name:	c-Fgr Antibody
Catalog Number:	B8105
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FGR. AA range:381-430
Concentration:	1 mg/ml
Gene Name:	FGR
Protein Name:	Tyrosine-protein kinase Fgr
Other Name:	FGR; SRC2; Tyrosine-protein kinase Fgr; Gardner-Rasheed feline sarcoma viral; v-fgr) oncogene homolog; Proto-oncogene c-Fgr; p55-Fgr; p58-Fgr; p58c-Fgr
Human Swiss Prot No.:	P09769
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	c-Fgr Polyclonal Antibody detects endogenous levels of c-Fgr protein.
Storage and Stability:	-20°C/1 year