

Antibodies Product Datasheet

Product Name:	Mnk1 Antibody
Catalog Number:	B8143
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MNK1. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	MKNK1
Protein Name:	MAP kinase-interacting serine/threonine-protein kinase 1
Other Name:	MKNK1; MNK1; MAP kinase-interacting serine/threonine-protein kinase 1; MAP kinase signal-integrating kinase 1; MAPK signal-integrating kinase 1; Mnk1
Human Swiss Prot No.:	Q9BUB5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mnk1 Polyclonal Antibody detects endogenous levels of Mnk1 protein.
Storage and Stability:	-20°C/1 year