

Antibodies Product Datasheet

Product Name:	SIK1 Antibody
Catalog Number:	B8185
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SIK. AA range:148-197
Concentration:	1 mg/ml
Gene Name:	SIK1
Protein Name:	Serine/threonine-protein kinase SIK1
Other Name:	SIK1; SIK; SNF1LK; Serine/threonine-protein kinase SIK1; Salt-inducible kinase 1; SIK-1; Serine/threonine-protein kinase SNF1-like kinase 1; Serine/threonine-protein kinase SNF1LK
Human Swiss Prot No.:	P57059
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SIK1 Polyclonal Antibody detects endogenous levels of SIK1 protein.
Storage and Stability:	-20°C/1 year