

## Antibodies Product Datasheet

|                               |                                                                                                                                                            |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | MLK3 Antibody                                                                                                                                              |
| <b>Catalog Number:</b>        | B8375                                                                                                                                                      |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                          |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA                                                                                                                         |
| <b>Dilutions:</b>             | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.                                                                 |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                           |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                          |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                      |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human MLK3. AA range:640-689                                                           |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                    |
| <b>Gene Name:</b>             | MAP3K11                                                                                                                                                    |
| <b>Protein Name:</b>          | Mitogen-activated protein kinase kinase kinase 11                                                                                                          |
| <b>Other Name:</b>            | MAP3K11; MLK3; PTK1; SPRK; Mitogen-activated protein kinase kinase kinase 11; Mixed lineage kinase 3; Src-homology 3 domain-containing proline-rich kinase |
| <b>Human Swiss Prot No.:</b>  | Q16584                                                                                                                                                     |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                    |
| <b>Specificity:</b>           | MLK3 Polyclonal Antibody detects endogenous levels of MLK3 protein.                                                                                        |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                               |