

Antibodies Product Datasheet

Product Name:	PDX-1 Antibody
Catalog Number:	B8398
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PDX1. AA range:27-76
Concentration:	1 mg/ml
Gene Name:	PDX1
Protein Name:	Pancreas/duodenum homeobox protein 1
Other Name:	PDX1; IPF1; Pancreas/duodenum homeobox protein 1; PDX-1; Glucose-sensitive factor; GSF; Insulin promoter factor 1; IPF-1; Insulin upstream factor 1; IUF-1; Islet/duodenum homeobox-1; IDX-1; Somatostatin-transactivating factor 1; STF-1
Human Swiss Prot No.:	P52945
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PDX-1 Polyclonal Antibody detects endogenous levels of PDX-1 protein.
Storage and Stability:	-20°C/1 year