

Antibodies Product Datasheet

Product Name:	Lck BP-1 Antibody
Catalog Number:	B8462
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Lck BP-1 . at AA range: 340-420
Concentration:	1 mg/ml
Gene Name:	HCLS1
Protein Name:	Hematopoietic lineage cell-specific protein
Other Name:	HCLS1; HS1; Hematopoietic lineage cell-specific protein; Hematopoietic cell-specific LYN substrate 1; LckBP1; p75
Human Swiss Prot No.:	P14317
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Lck BP-1 Polyclonal Antibody detects endogenous levels of Lck BP-1 protein.
Storage and Stability:	-20°C/1 year