

Antibodies Product Datasheet

Product Name:	GABAB R2 Antibody
Catalog Number:	B8506
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from GABAB R2 . at AA range: 830-910
Concentration:	1 mg/ml
Gene Name:	GABBR2
Protein Name:	Gamma-aminobutyric acid type B receptor subunit 2
Other Name:	GABBR2; GPR51; GPRC3B; Gamma-aminobutyric acid type B receptor subunit 2; GABA-B receptor 2; GABA-B-R2; GABA-BR2; GABABR2; Gb2; G-protein coupled receptor 51; HG20
Human Swiss Prot No.:	O75899
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABAB R2 Polyclonal Antibody detects endogenous levels of GABAB R2 protein.
Storage and Stability:	-20°C/1 year