

Antibodies Product Datasheet

Product Name:	LAT Antibody
Catalog Number:	B9043
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from LAT . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	LAT
Protein Name:	Linker for activation of T-cells family member 1
Other Name:	LAT; Linker for activation of T-cells family member 1; 36 kDa phospho-tyrosine adapter protein; pp36; p36-38
Human Swiss Prot No.:	O43561
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LAT Polyclonal Antibody detects endogenous levels of LAT protein.
Storage and Stability:	-20°C/1 year