

Antibodies Product Datasheet

Product Name:	AP-2 α / β Antibody
Catalog Number:	C0127
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AP-2. AA range:388-437
Concentration:	1 mg/ml
Gene Name:	TFAP2A/TFAP2B
Protein Name:	Transcription factor AP-2-alpha/beta
Other Name:	TFAP2A; AP2TF; TFAP2; Transcription factor AP-2-alpha; AP2-alpha; AP-2 transcription factor; Activating enhancer-binding protein 2-alpha; Activator protein 2; AP-2; TFAP2B; Transcription factor AP-2-beta; AP2-beta; Activating enhancer-bind
Human Swiss Prot No.:	P05549/Q92481
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AP-2 α / β Polyclonal Antibody detects endogenous levels of AP-2 α / β protein.
Storage and Stability:	-20°C/1 year