

Antibodies Product Datasheet

Product Name:	Dysferlin Antibody
Catalog Number:	C0176
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Dysferlin. AA range:1981-2030
Concentration:	1 mg/ml
Gene Name:	DYSF
Protein Name:	Dysferlin
Other Name:	DYSF; FER1L1; Dysferlin; Dystrophy-associated fer-1-like protein; Fer-1-like protein 1
Human Swiss Prot No.:	O75923
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dysferlin Polyclonal Antibody detects endogenous levels of Dysferlin protein.
Storage and Stability:	-20°C/1 year