

Antibodies Product Datasheet

Product Name:	Factor VIII Antibody
Catalog Number:	C0187
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Factor VIII. AA range:2161-2210
Concentration:	1 mg/ml
Gene Name:	F8
Protein Name:	Coagulation factor VIII
Other Name:	F8; F8C; Coagulation factor VIII; Antihemophilic factor; AHF; Procoagulant component
Human Swiss Prot No.:	P00451
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Factor VIII Polyclonal Antibody detects endogenous levels of Factor VIII protein.
Storage and Stability:	-20°C/1 year