

Antibodies Product Datasheet

Product Name:	GFR α -1 Antibody
Catalog Number:	C0200
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GFR α -1. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	GFRA1
Protein Name:	GNDF family receptor α -1
Other Name:	GFRA1; GDNFRA; RETL1; TRNR1; GDNF family receptor α -1; GDNF receptor α -1; GDNFR- α -1; GFR- α -1; RET ligand 1; TGF- β -related neurotrophic factor receptor 1
Human Swiss Prot No.:	P56159
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GFR α -1 Polyclonal Antibody detects endogenous levels of GFR α -1 protein.
Storage and Stability:	-20°C/1 year