

Antibodies Product Datasheet

Product Name:	Granzyme B/H Antibody
Catalog Number:	C0215
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Granzyme B. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	GZMB/GZMH
Protein Name:	Granzyme B/H
Other Name:	GZMB; CGL1; CSPB; CTLA1; GRB; Granzyme B; C11; CTLA-1; Cathepsin G-like 1; CTSG1; Cytotoxic T-lymphocyte proteinase 2; Lymphocyte protease; Fragmentin-2; Granzyme-2; Human lymphocyte protein; HLP; SECT; T-cell serine protease 1-3E; GZMH; C
Human Swiss Prot No.:	P10144/P20718
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Granzyme B/H Polyclonal Antibody detects endogenous levels of Granzyme B/H protein.
Storage and Stability:	-20°C/1 year