

Antibodies Product Datasheet

Product Name:	Integrin β 5 Antibody
Catalog Number:	C0235
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Integrin beta5. AA range:721-770
Concentration:	1 mg/ml
Gene Name:	ITGB5
Protein Name:	Integrin beta-5
Other Name:	ITGB5; Integrin beta-5
Human Swiss Prot No.:	P18084
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin β 5 Polyclonal Antibody detects endogenous levels of Integrin β 5 protein.
Storage and Stability:	-20°C/1 year