

Antibodies Product Datasheet

Product Name:	PIASx Antibody
Catalog Number:	C0362
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PIAS2. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	PIAS2
Protein Name:	E3 SUMO-protein ligase PIAS2
Other Name:	PIAS2; PIASX; E3 SUMO-protein ligase PIAS2; Androgen receptor-interacting protein 3; ARIP3; DAB2-interacting protein; DIP; Msx-interacting zinc finger protein; Miz1; PIAS-NY protein; Protein inhibitor of activated STAT x; Protein inhibitor
Human Swiss Prot No.:	O75928
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PIASx Polyclonal Antibody detects endogenous levels of PIASx protein.
Storage and Stability:	-20°C/1 year