

Antibodies Product Datasheet

Product Name:	TNF-R1 Antibody
Catalog Number:	C0375
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNF Receptor I. AA range:381-430
Concentration:	1 mg/ml
Gene Name:	TNFRSF1A
Protein Name:	Tumor necrosis factor receptor superfamily member 1A
Other Name:	TNFRSF1A; TNFAR; TNFR1; Tumor necrosis factor receptor superfamily member 1A; Tumor necrosis factor receptor 1; TNF-R1; Tumor necrosis factor receptor type I; TNF-RI; TNFR-I; p55; p60; CD antigen CD120a
Human Swiss Prot No.:	P19438
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TNF-R1 Polyclonal Antibody detects endogenous levels of TNF-R1 protein.
Storage and Stability:	-20°C/1 year