

Antibodies Product Datasheet

Product Name:	Anoctamin-9 Antibody
Catalog Number:	C10001
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TM16J. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	ANO9
Protein Name:	Anoctamin-9
Other Name:	ANO9; PIG5; TMEM16J; TP53I5; Anoctamin-9; Transmembrane protein 16J; Tumor protein p53-inducible protein 5; p53-induced gene 5 protein A1A5B4
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Anoctamin-9 Polyclonal Antibody detects endogenous levels of Anoctamin-9 protein.
Storage and Stability:	-20°C/1 year